

Comparison of dNTPs Mix from Vivantis and Supplier P

Amplification of MCS region with inserts 1.4kb from pTZ using M13 primer in a 50 μ l reaction.

5 μ l of PCR product is loaded per lane and electrophoresed in 1.0% TBE agarose gel.

2 different brands of dNTPs mixes are tested in 0.2mM, 0.1mM and 0.05mM in a 50 μ l PCR reaction.

Experimental Date: 25th July 2016




Figure 1: 10mM dNTPs mix is diluted into 2mM, 1mM and 0.5mM. Amplification of pTZ template using M13 primer and 5 μ l of different concentration of dNTPs mix. 5 μ l of PCR product is loaded per lane. Expected PCR product size amplified is 1.4kb.

Legend:

- 1: Amplification using 0.3ng/ μ l of pTZ DNA template, 0.2 μ M of primers and 0.2mM dNTPs mix
- 2: Amplification using 0.3ng/ μ l of pTZ DNA template, 0.2 μ M of primers and 0.1mM dNTPs mix
- 3: Amplification using 0.3ng/ μ l of pTZ DNA template, 0.2 μ M of primers and 0.05mM dNTPs mix
- ve: Amplification using 0.2 μ M of primers and 0.2mM dNTPs mix
- M: VC 1kb DNA ladder



Conclusion:

Both dNTPs mixes are performed up to standard; the performance of Vivantis's dNTPs mix is equivalent to the performance of Supplier P's dNTPs mix.

Prepared by,
Vivantis Technical Team
25th July 2016

Pairing Nature with Scientific Discoveries

Vivantis Technologies Sdn Bhd 587389-D

Revongen Corporation Center

No12A, Jalan TP5, Taman Perindustrian UEP,
47600 Subang Jaya, Selangor Darul Ehsan, Malaysia.

E : info@vivantechnologies.com

W: www.vivantechnologies.com

T: +6 03 8025 1603

F: +6 03 8025 1637/1354