



A typical freezer in a molecular biology laboratory consists of a host of the quintessential restriction enzymes from various sources. Next come the laborious hunt for different buffer charts to culminate in a compromise of buffers for double digestions. Vivantis presents a solution to this task! Our 10-go kits consist of selected enzymes for your use, coupled with recommended optimal buffers for double digestions. This is complemented with an unbelievable price performance ratio with the offer of "buy 1, get more". Our 10-go kits present sets of restriction enzymes to provide you with hassle free tools for your molecular biology applications, saving both time and money!

10-go, giving scientists time to tango!

### 10-GO RE KIT SET (C) - RK3000

Kit Component: Restriction Endonucleases (10 types)

1ml Diluent Viva Buffer A

1ml 10X Buffers (Buffer V1, V2, V3, V4, V5, Universal Buffer, and Buffer *EcoR* I, Buffer *Bsp19* I)

(BSA included in all Reaction Buffer)

LOT #

STORE AT -20°C

### RESTRICTION ENDONUCLEASES SPECIFICATION:

ENZYME	(u/μL)	TOTAL UNIT	REG.SITE	OPT. BUFFER	OPT.(°C)	THERMA INAC.	ACTIVITY IN REACTION BUFFER				
							V1	V2	V3	V4	V5
<i>Apa</i> I	20	150	5'..GGGCC^C..3'	V5	37	65°C, 20min.	75%	75%	75%	75%	100%
<i>Bam</i> H I	20	400	5'..G^GATCC..3'	UB1X	37	65°C, 20min.	75%	75%	50%	75%	50%
<i>Bsp19</i> I ( <i>Nco</i> I)	20	100	5'..C^CATGG..3'	SP	37	65°C, 20min.	25%	50%	75%	75%	50%
<i>Cci</i> N I ( <i>Not</i> I)	5	50	5'..GC^GGCCGC..3'	V5	37	65°C, 20min.	50%	75%	75%	75%	100%
<i>Eco</i> R I	20	500	5'..G^AATTC..3'	SP	37	65°C, 20min.	50%	50%	100%	100%	50%
<i>Hind</i> III	20	400	5'..A^AGCTT..3'	V2	37	65°C, 20min.	75%	100%	75%	75%	75%
<i>Kpn</i> I	20	200	5'..GGTAC^C..3'	V1	37	Non-Thermal Inac.	100%	25%	25%	25%	75%
<i>Psp124B</i> I ( <i>Sac</i> I)	20	100	5'..GAGCT^C..3'	V2	37	65°C, 20min.	100%	100%	75%	100%	100%
<i>Sfr303</i> I ( <i>Sac</i> II)	5	100	5'..CCGC^GG..3'	V1	37	65°C, 20min.	100%	75%	75%	50%	50%
<i>Xba</i> I	20	300	5'..T^CTAGA..3'	V5	37	65°C, 20min.	10%	75%	75%	10%	100%

#### REMARK:

*Apa* I Blocked by overlapping dcm-methylation (CC<sup>m</sup>WGG): GGGCC<sup>m</sup>WGG  
*Bam*H I High enzyme concentration may result in Star Activity.  
*Cci*N I High enzyme concentration may result in Star Activity.  
*Eco*R I High enzyme concentration may result in Star Activity.  
*Xba* I Blocked by overlapping dam-methylation (G<sup>m</sup>ATC): TCTAG<sup>m</sup>ATC

#### BUFFER COMPOSITION:

<b>V1<sub>eff</sub></b>	10mM Tris-HCl ( pH 7.5 at 30°C), 10mM MgCl <sub>2</sub> , and 100μg/ml BSA.
<b>V2<sub>eff</sub></b>	10mM Tris-HCl ( pH 7.5 at 30°C), 10mM MgCl <sub>2</sub> , 50mM NaCl, and 100μg/ml BSA.
<b>V3<sub>eff</sub></b>	50mM Tris-HCl ( pH 7.5 at 30°C), 10mM MgCl <sub>2</sub> , 100mM NaCl, and 100μg/ml BSA.
<b>V4<sub>eff</sub></b>	10mM Tris-HCl ( pH 8.5 at 30°C), 10mM MgCl <sub>2</sub> , 100mM KCl, and 100μg/ml BSA.
<b>V5<sub>eff</sub></b>	30mM Tris-acetate ( pH 7.9 at 30°C), 10mM Mg-acetate, 60mM K-acetate, and 100μg/ml BSA.
<b>UB<sub>eff</sub></b>	1.0X UB 25mM Tris-acetate (pH 7.6 at 30°C), 10mM Mg-acetate, 100mM K-acetate, 7mM 2-mercaptoethanol and 50μg/ml BSA
<b>SP<sub>eff</sub></b>	Buffer <i>EcoR</i> I 500mM Tris-HCl (pH 7.5 at 30°C), 10mM MgCl <sub>2</sub> , 100mM NaCl, 0.02% triton X-100, and 0.1mg / ml BSA.

### BUFFER COMPATIBILITY CHART

RK3000	<i>Apa</i> I	<i>Bam</i> H I	<i>Bsp19</i> I	<i>Cci</i> N I	<i>Eco</i> R I	<i>Hind</i> III	<i>Kpn</i> I	<i>Psp124B</i> I	<i>Sfr303</i> I	<i>Xba</i> I
<i>Apa</i> I		V2/V5	V5	V5	V3/V4	V2/V5	V1/V5	V5	V1	V5
<i>Bam</i> H I	V2/V5		V2-V5	V2/V5	V4	V1/V4/V5	V1	V1/V4/V5	V1/V2	V2/V5
<i>Bsp19</i> I	V5	V2-V5		V5	V3/V4	V2-V5	V1/V5	V4/V5	V3	V5
<i>Cci</i> N I	V5	V2/V5	V5		V3/V4	V2/V5	V5	V5	V2/V3	V5
<i>Eco</i> R I	V3/V4	V4	V3/V4	V3/V4		V3/V4	V1	V4	V3	1.5X UB
<i>Hind</i> III	V2/V5	V1/V4/V5	V2-V5	V2/V5	V3/V4		V1	V2	V1/V2	V2/V5
<i>Kpn</i> I	V1/V5	V1	V1/V5	V5	V1	V1		V1	V1	V5
<i>Psp124B</i> I	V5	V1/V4/V5	V4/V5	V5	V4	V2	V1		V1	V5
<i>Sfr303</i> I	V1	V1/V2	V3	V2/V3	V3	V1/V2	V1	V1		V2/V3
<i>Xba</i> I	V5	V2/V5	V5	V5	1.5X UB	V2/V5	V5	V5	V2/V3	

#### Activity:

	100%		Specific buffer
	75%		Not applicable
	50%		

ENZYME	LIGATION / RECURTING ASSAY	OVERDIGESTION ASSAY
<b>Apa I</b>	After 20-fold overdigestion with <b>Apa I</b> , more than 90% of the DNA fragments can be ligated and recut.	An unaltered banding pattern was observed after 1µg of DNA was digested with 40u of <b>Apa I</b> for 16 hours at 37°C (Without BSA).
<b>BamH I</b>	After 20 -fold overdigestion with <b>BamH I</b> , more than 90% of the DNA fragments can be ligated and recut.	An unaltered banding pattern was observed after 1µg of DNA was digested with 40 u of <b>BamH I</b> for 16 hours at 37°C .
<b>Bsp19 I</b> ( <i>Nco I</i> )	After 20-fold overdigestion with <b>Bsp19 I</b> , more than 90% of DNA fragments can be ligated and recut.	An unaltered banding pattern was observed after 1µg of DNA was digested with 40u of <b>Bsp19 I</b> for 16 hours at 37°C (Without BSA).
<b>CciN I</b> ( <i>Not I</i> )	After 5-fold overdigestion with <b>CciN I</b> , more than 90% of the DNA fragments can be ligated and recut.	An unaltered banding pattern was observed after 1µg of DNA was digested with 10 u of <b>CciN I</b> for 16 hours at 37°C (Without BSA).
<b>EcoR I</b>	After 20-fold overdigestion with <b>EcoR I</b> , more than 95% of the DNA fragments can be ligated and recut.	An unaltered banding pattern was observed after 1µg of DNA was digested with 40u of <b>EcoR I</b> for 16 hours at 37°C (Without BSA).
<b>Hind III</b>	After 20-fold overdigestion with <b>Hind III</b> , more than 90% of the DNA fragments can be ligated and recut.	An unaltered banding pattern was observed after 1µg of DNA was digested with 40u of <b>Hind III</b> for 16 hours at 37°C.
<b>Kpn I</b>	After 20-fold overdigestion with <b>Kpn I</b> , more than 90% of the DNA fragments can be ligated and recut.	An unaltered banding pattern was observed after 1µg of DNA was digested with 40u of <b>Kpn I</b> for 16 hours at 37°C (Without BSA).
<b>Psp124B I</b>	After 20-fold overdigestion with <b>Psp124B I</b> , more than 95% of the DNA fragments can be ligated and recut.	An unaltered banding pattern was observed after 1µg of DNA was digested with 40u of <b>Psp124B I</b> for 16 hours at 37°C (Without BSA).
<b>Sfr303 I</b>	After 5-fold overdigestion with <b>Sfr303 I</b> , more than 90% of the DNA fragments can be ligated and recut.	An unaltered banding pattern was observed after 1µg of DNA was digested with 10u of <b>Sfr303 I</b> for 16 hours at 37°C.
<b>Xba I</b>	After 20-fold overdigestion with <b>Xba I</b> , 90% of the DNA fragments can be ligated and recut.	An unaltered banding pattern was observed after 1µg of DNA was digested with 40u of <b>Xba I</b> for 16 hours at 37°C.