

1<sup>st</sup> July 2018

Dear valued customers,

**CHANGE OF PRODUCT SPECIFICATION NOTICE**  
**GF-1 Soil Sample DNA Extraction Kit**

This letter serves as a formal notification that Vivantis Technologies Sdn Bhd has decided to change specification of **GF-1 Soil Sample DNA Extraction Kit** with effective from **1<sup>st</sup> July 2018**.

<b>Changes</b>	<b>Current Specification</b>	<b>New Specification</b>
Glass beads	0.5mm	Mixture of 0.5mm and 0.1mm
Buffers	HTR Reagent Buffer DS Buffer SLX Mlus P2 Buffer XP1 Buffer SPW Wash Buffer Elution Buffer	Buffer SL1 Buffer SL2 Buffer SB SPW Wash Buffer 1 SPW Wash Buffer 2 Elution Buffer
Columns	Blue-color columns	White-color columns

**Ordering Information:**

<b>Catalogue No.</b>	<b>Description</b>	<b>Pack Size</b>
GF-SD-005	GF-1 Soil Sample DNA Extraction Kit	5 preps
GF-SD-025	GF-1 Soil Sample DNA Extraction Kit	25 preps

Date of termination of current specification of GF-1 Soil Sample DNA Extraction Kit:  
 As soon as the inventory held by the Company depleted or estimated by **31<sup>st</sup> July 2018**.

Date of ordering of new specification of GF-1 Soil Sample DNA Extraction Kit: **16<sup>th</sup> July 2018**.

**Last time purchase conditions:**

1. Date for last order acceptance for old spec of GF-1 Soil Sample DNA Extraction Kit: **31<sup>st</sup> July 2018**.
2. The order shall not be fulfilled if there is no inventory held by the Company.
3. All purchase orders are non-cancellable and non-returnable.

**Pairing Nature with Scientific Discoveries**



We sincerely apologize for any inconvenience that may have caused to our customers. As regards this decision made by the Company, please note that there are no quality or safety issues with the existing product. Your kind understanding in this matter is greatly appreciated.

To provide the best product performance for your research application, the Company have supplied of the improved version of GF-1 Soil Sample DNA Extraction Kit since July 2018 as an alternative and replacement for current specification of product. With support and positive feedback from customers, we are confident of the performance of GF-1 Soil Sample DNA Extraction Kit with new specification.

We sincerely value and appreciate the business relationship established with your esteemed organization and we look forward to your continuous interest and support in our products. We fervently anticipate to fulfilling the requirement of upcoming new products with high standard of quality.

For further questions and information, we are your disposal at any time. Please contact our sales representatives on outstanding issues at +603 8025 1603. You can also reach us at [info@vivantechnologies.com](mailto:info@vivantechnologies.com).

We appreciate your business and look forward to assisting you.

Yours sincerely,

Vivantis Technologies Sdn Bhd

**Pairing Nature with Scientific Discoveries**

Vivantis Technologies Sdn Bhd 587389-D

Revongen Corporation Center

No12A, Jalan TP5, Taman Perindustrian UEP,  
47600 Subang Jaya, Selangor Darul Ehsan, Malaysia.

E : [info@vivantechnologies.com](mailto:info@vivantechnologies.com)  
W: [www.vivantechnologies.com](http://www.vivantechnologies.com)

T: +6 03 8025 1603  
F: +6 03 8025 1637/1354