

## pBR322 DNA

Product No : NN1404  
Quantity : 100µg

Lot :  
Expiry date :  
Concentration : 0.5µg/µl

Store at -20°C [info@vivantechnologies.com](mailto:info@vivantechnologies.com)

### Description:

pBR322 DNA is a commonly used plasmid cloning vector. Isolated from *E. coli* strain Top10F<sup>1</sup>.

### Storage Buffer:

10mM Tris-HCl (pH 8.0) and 1mM EDTA.

### Quality Control:

Gel analysis for purity. Digested with *Alu I*, *Bsn I* and *Hind III* for cut pattern confirmation.



Product Use Limitation

This product is for research purposes and in vitro use only.

v i v a n t i s | [www.vivantechnologies.com](http://www.vivantechnologies.com)

DSNN1404\_rev0\_010311