

pUC19 DNA

Product No : NN1406
Quantity : 50µg

Lot :
Expiry date :
Concentration : 0.5µg/µl

Store at -20°C

info@vivantechnologies.com

Description:

pUC19 is commonly used plasmid cloning vectors. Isolated from *E. coli* strain Top10F⁺.

Storage Buffer:

10mM Tris-HCl (pH 8.0) and 1mM EDTA.

Quality Control:

Gel analysis for purity. Digested with **EcoR I**, **Hind III** and **Msp I** for cut pattern confirmation.



Product Use Limitation

This product is for research purposes and *in vitro* use only.